Medtronic

FREESTYLE®

Aortic Root Bioprosthesis / Bioprothèse de racine aortique / Aortenwurzelbioprothese / Bioprótesis de la raíz aórtica / Aortawortelbioprothese / Bioprotesi della radice aortica / Aortarotbioproteze / Aortan tyven bioproteesi / Biologisk aortarotklaff / Βιοπρόσθεση αορτικής ρίζας / Aortarod-bioproteze / Biologiczna zastawka aortalna wraz z opuszką aorty / Bioprótese de raiz aórtica / Aortik Kõk Bioprotezi / Bioprotéza aortálního kořene / Aortagyök biológiai műbillentyű / Bioprotéza aortálneho koreňa / Aordi lähte bioprotees / Биопротеза за аортния корен / Bioproteza korijena aorte / Aortas saknes bioprotēze / Aortos šaknies bioprotezas / Bioproteză pentru rădăcină aortică / Биопротез корня аорты / Bioproteza korena aorte / Bioproteza za koren aorte

Instructions for Use / Mode d'emploi / Gebrauchsanweisung / Instrucciones de uso / Gebruiksaanwijzing / Istruzioni per l'uso / Bruksanvisning / Käyttöohjeet / Bruksanvisning / Οδηγίες χρήσης / Brugsanvisning / Instrukcja użytkowania / Instruções de utilização / Kullanım Talimatları / Návod k použití / Használati útmutató / Pokyny na používanie / Kasutusjuhend / Инструкции за употреба / Upute za upotrebu / Lietošanas pamācība / Naudojimo instrukcijos / Instrucțiuni de utilizare / Инструкция по эксплуатации / Uputstva za upotrebu / Navodila za uporabo

Explanation of symbols on package labeling / Explication des symboles des étiquettes sur l'emballage / Erläuterung der Symbole auf dem Verpackungsetikett / Explicación de los símbolos en la documentación del envase / Verklaring van de symbolen op de verpakkingslabels / Spiegazione dei simboli sulle etichette della confezione / Forklaring av symboler på emballasjen / Pakkauksen symbolien selitykset / Förklaring av symboler på förpackningsetiketten / Επεξήγηση των συμβόλων στη σήμανση της συσκευασίας / Forklaring af symboler på emballagens mærkater / objaśnienie symboli znajdujących się na etykietach opakowania / Explicação dos símbolos na documentação da embalagem / Ambalaj etiketi üzerindeki sembollerin açıklaması / Vysvětlení symbolů uvedených na štítcích obalu / A csomagoláson látható szimbólumok jelentése / Vysvetlivky k symbolom na označení balenia / Pakendi siltidel olevate tingmärkide tähendus / Объяснение на символите от етикетите на опаковката / Objašnjenje simbola na etiketi na ambalaži / Uz iepakojuma esošo simbolu izskaidrojums / Pakuotės etiketės nurodytų simbolių reikšmės / Explicarea simbolurilor cuprinse în etichetele de pe ambalaj / Объяснение символов, приведенных на упаковке / Objašnjenje simbola na oznakama na pakovanju / Razlaga simbolov na embalaži

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
1 DEVICE DESCRIPTION

The FREESTYLE Aortic Root Bioprosthesis (FREESTYLE Bioprosthesis), Model 995, consists of a porcine aortic root preserved in buffered 0.2% glutaraldehyde with a thin polyester cloth covering added to strengthen the proximal (inflow) suture line and to cover any exposed porcine myocardium.

Surgeon's flags, located 120 degrees apart, are provided at the inflow aspect of the bioprosthesis to facilitate uniform placement of stitches for the proximal suture line. Colored stitching, circumferentially placed around the cloth cover, demarcates the upper limit for proximal implantation sutures. The size of the bioprosthesis is the outside diameter at the inflow edge. The FREESTYLE Aortic Root Bioprostheses are available in the sizes shown in Table 1. Profile and inner cloth height dimensions are also provided in the table.

Table 1. FREESTYLE Full Root Aortic Valve, Model 995
Available Sizes and Dimensions

Size (mm)	A Outside Diameter (+0.5 mm -0.0 mm)	B Profile Height (±2 mm)	C Inner Cloth Height (±0.5 mm)
19	19.0	30	3.0
21	21.0	32	3.0
23	23.0	32	3.0
25	25.0	34	3.0
27	27.0	34	3.0
29	29.0	39±3 mm	3.0

The FREESTYLE Bioprosthesis is treated with an alpha amino oleic acid antimineralization process, AOA, which has been shown to mitigate porcine leaflet calcification in animal studies.

This device contains no metals and, therefore, poses no known hazards in all MR environments.

2 INDICATIONS

The FREESTYLE Bioprostheses are indicated for the replacement of pathologic or prosthetic aortic valves.

3 CONTRAINDICATIONS

Patients with chronic renal failure, those undergoing hemodialysis, or those with abnormal calcium metabolism, may have an increased tendency toward calcification of valvular bioprostheses. Bioprosthetic valves are also susceptible to calcification in children. Patients with Marfan's syndrome or chronic endocarditis may be contraindicated for a bioprosthesis. Accordingly, the physician should carefully consider these potential hazards when selecting an appropriate valve substitute for such patients.

4 WARNINGS

This device was designed for single patient use only. Do not reuse, reprocess, or resterilize this product. Reuse, reprocessing, or resterilization may compromise the structural integrity of the device and/or create a risk of contamination of the device, which could result in patient injury, illness, or death.

Check the shipping temperature indicator inside the carton. If the shipping temperature indicator window is black, the valve is not suitable for clinical use.

DO NOT RESTERILIZE THE VALVE BY ANY METHOD. Exposure of the bioprosthesis and container to irradiation, steam, ethylene oxide or other chemical sterilants will render the bioprosthesis completely unfit for use.

Do not use any valve that has been damaged. Extreme care must be taken to prevent damage to the delicate valve tissue. Do not handle the leaflet portion of the bioprosthesis with instruments. Even the most minor perforation may enlarge in time to produce significant impairment of valve function. Should a valve be damaged during insertion, do not attempt repair.

Do not rinse the bioprosthesis in antibiotic solution.

5 PRECAUTIONS

To prevent drying of the valve tissue, periodically irrigate the bioprosthesis with sterile normal saline during implantation.

Identify the inflow and outflow of the bioprosthesis prior to suturing.

When selecting a valve size, consideration of the cardiac anatomy is necessary, and care must also be taken to select a valve that adequately provides for the hemodynamic requirements of the patient.

Do not use cutting needles, as they may cause structural damage to the bioprosthesis. The potential for damage to the bioprosthesis should be considered before passing catheters, surgical instruments or transvenous pacing leads across the valve.

Passage of a catheter through any bioprosthesis may damage the valve and is, therefore, not recommended.

Though the coronary arteries of the FREESTYLE Bioprosthesis are ligated during the manufacturing process, leakage from the area of the coronary ligation may still occur. In any procedure for which one of the coronary arteries are retained, reinforcement of the ligation may be required to prevent leakage.

Suture ends must be trimmed close to the knot to prevent abrasion of leaflet tissue.

Patients exhibiting an absence of sinus rhythm, a hypertrophied left atrium, or having a history of previous thromboembolic complications should be considered for indefinite anticoagulation therapy.

Caution: Do not invert the bioprosthesis when suturing. Inversion may result in elongated suture holes, tears, and/or distortion leading to stenosis and incompetence.

DO NOT UNDERSIZE THE FREESTYLE BIOPROSTHESIS. Unlike fresh, cryopreserved, or antibiotic sterilized aortic homografts, the leaflets of the FREESTYLE Bioprosthesis may not stretch to adequately compensate for underestimation of supracommissural diameter during diastole. Thus if the indicated annular/supracommissural size falls between consecutive implantation diameters, the larger size bioprosthesis should be implanted.

Note: Care should be taken to ensure that excessive mismatch of the annular diameter and the supracommissural ridge diameter does not occur. Excessive mismatch may cause central incompetence and/or stenosis of the bioprosthesis.

Use extreme caution if tailoring the bioprosthesis to fit the anatomical requirements of a particular patient's coronary sinuses or ostia. Improper trimming may result in immediate or delayed damage to and/or dysfunction of the bioprosthesis.

6 POTENTIAL ADVERSE EVENTS

Adverse events potentially associated with the use of bioprosthetic heart valves include:

- cardiac dysrhythmias
- death
- endocarditis
- hemolysis
- hemorrhage, anticoagulant/antiplatelet-related
- leak, transvalvular or paravalvular
- nonstructural dysfunction (pannus, suture, inappropriate sizing)
- structural deterioration (calcification, leaflet tear, intracuspal hematoma, pseudo aneurysm, or other)
- thromboembolism
- valve thrombosis
- root dilatation

7 INDIVIDUALIZATION OF TREATMENT

Long-term anticoagulant and/or antiplatelet therapy should be considered in patients with a dilated left atrium, a history of thromboembolic events, or a cardiac rhythm of atrial fibrillation or atrial flutter.

8 PATIENT COUNSELING INFORMATION

Patients may require anticoagulation and/or antiplatelet therapy for an indefinite period based on each patient's condition. Patients with bioprostheses are at risk for bacteremia (eg, undergoing dental procedures) and should be advised about prophylactic antibiotic therapy.

Patients should be encouraged to carry the Implanted Device Identification Card, provided by Medtronic, with them at all times.

9 HOW SUPPLIED

9.1 Packaging

The FREESTYLE Bioprosthesis is chemically sterilized and is supplied STERILE in a buffered 0.2% glutaraldehyde storage solution. Sterility is compromised if the glass jar and lid container is opened and/or damaged. The outside of the container is NOT sterile and should not be placed in the sterile field.

9.2 Storage

The FREESTYLE Bioprosthesis must be stored between 5°C and 25°C (41°F and 77°F). Refrigeration is not required, and freezing may damage the bioprosthesis. Room temperature storage up to 25°C (77°F) is satisfactory provided the bioprosthesis is not exposed to sunlight or other ultraviolet light sources or placed where significant temperature fluctuations may occur.

Appropriate inventory control should be maintained so that bioprostheses with earlier use by dates are preferentially implanted and expiration is avoided.

10 INSTRUCTIONS FOR USE

10.1 Handling and Preparation Instructions

Proper bioprosthesis size selection is a very important part of heart valve replacement. The internal diameter of the patient's aortic root at the annulus and supracommissural areas may be measured preoperatively, during diastole, using angiographic and/or echocardiographic techniques. The size selection of a FREESTYLE Bioprosthesis is aided by the use of FREESTYLE Aortic Obturators (Model 7990). Use only FREESTYLE Aortic Obturators to select the appropriate size bioprosthesis. For further information, refer to the FREESTYLE Aortic Obturators Instructions for Use.

Within the sterile operative field, prepare 3 rinse basins, each containing 500 mL of sterile normal saline solution. The exterior of the device container and lid are NONSTERILE. Examine the tamper strips to verify that the container has not been damaged or previously opened. Do not use if *all* the tamper strips are damaged. Turn the lid counterclockwise and open the container (Figure 1).

The bioprosthesis and all internal packaging components within the container are STERILE and must be handled accordingly. With the thumb and index finger, grasp the retainer and slowly lift it out of the container allowing for drainage of the glutaraldehyde storage solution (Figure 2).

Snap open the cap from the retainer body and place the valve directly into the free hand (Figures 3 and 4).

Have an assistant record the valve identification number of the bioprosthesis in the patient's record.

Carefully cut the identification tag from the valve and discard the tag (Figure 5).

Note: Be careful not to cut the cloth or valve tissue when removing the identification tag. Remove any remnants of the identification tag suture from the valve.

Submerge the valve into the first rinse basin. Do not touch the valve leaflets or squeeze the valve during the rinsing procedure. Gently swirl the valve in the solution for a minimum of 2 minutes in each of the 3 previously prepared rinse basins (Figure 6). The valve should remain in the third rinse basin until required by the surgeon.

10.2 Device Implantation

The FREESTYLE Bioprosthesis is designed for freehand implantation.

Care should be exercised when placing sutures through the sewing rim and aortic wall to prevent stitching through or perforation of the valve cusps. The colored stitches at the inflow demarcate the limited area for placing sutures in the sewing rim. Sutures should only be placed proximal to this demarcation line.

During implantation, periodically irrigate the valve with sterile normal saline to prevent drying of delicate valve tissue.

Do not use cutting needles as they could cause structural damage to the bioprosthesis.

Take care not to evert (roll outward) the inflow end of the bioprosthesis when suturing the valve to the patient's annulus. Eversion could damage the valve tissue.

Root Inclusion or Subcoronary Technique

In some patients, the height of the cloth covering the right muscle bar may exceed the height of the origin of the patient's right coronary. Forcing the bioprosthesis to fit in this situation could cause buckling or distortion of the porcine leaflets and/or aortic wall. Rotation of the bioprosthesis in such cases may be necessary. In the event of rotation, the noncoronary sinus must be trimmed. Use caution when trimming the noncoronary sinus to avoid damage to the valve leaflets.

If the bioprosthetic root is tailored to accommodate the patient's ostia, commissure posts should be fixed to the patient's aortic root 2 to 3 mm above the native commissure attachments to prevent leaflet prolapse and aortic regurgitation.

10.3 Accessories

Use only FREESTYLE Aortic Obturators (Model 7990) and the Medtronic Valve Handle (Model 7639) to determine the appropriate FREESTYLE Bioprosthesis size.

Caution: Accessories should be sterilized according to the Instructions for Use prior to use.

Caution: Do not use other manufacturers' valve obturators, or obturators for another Medtronic prosthesis to size the FREESTYLE Bioprosthesis.

10.4 Return of Explanted Bioprostheses

Medtronic is interested in obtaining recovered FREESTYLE Porcine Bioprostheses. When determined to be appropriate, explants will be studied by a consulting pathologist. A written report summarizing the findings will be returned to the physician. Product return kits, including an explant information form, are available by contacting Medtronic distribution centers or a Medtronic Sales Representative. It is important that the explant form be completely filled out. If a kit is not available, place the explanted bioprosthesis in a container of glutaraldehyde or 10% buffered formalin immediately after excision. For further instructions on the return of an explanted device, contact a Medtronic Sales Representative.

11 PATIENT INFORMATION

11.1 Registration Information

Note: Patient registration does not apply in countries where patient privacy laws conflict with providing patient information, including countries from the EU.

A Patient Registration Form is included in each device package. After implantation, please complete all requested information. The serial number may be found on the package and on the identification tag attached to the bioprosthesis. Return the original form to the Medtronic address indicated on the form and provide the temporary identification card to the patient prior to discharge.

An Implanted Device Identification Card is provided to the patient. The card contains the name and telephone number of the patient's physician as well as information that medical personnel would require in the event of an emergency.

12 DISCLAIMER OF WARRANTY

THE FOLLOWING DISCLAIMER OF WARRANTY APPLIES TO CUSTOMERS OUTSIDE THE UNITED STATES:

DISCLAIMER OF WARRANTY

ALTHOUGH THE FREESTYLE BIOPROSTHESIS, MODEL 995, HEREAFTER REFERRED TO AS “PRODUCT”, HAS BEEN CAREFULLY DESIGNED, MANUFACTURED AND TESTED PRIOR TO SALE, THE PRODUCT MAY FAIL TO PERFORM ITS INTENDED FUNCTION SATISFACTORILY FOR A VARIETY OF REASONS. THE WARNINGS CONTAINED IN THE PRODUCT LABELING PROVIDE MORE DETAILED INFORMATION AND ARE CONSIDERED AN INTEGRAL PART OF THIS DISCLAIMER OF WARRANTY. MEDTRONIC, THEREFORE, DISCLAIMS ALL WARRANTIES, BOTH EXPRESS AND IMPLIED, WITH RESPECT TO THE PRODUCT. MEDTRONIC SHALL NOT BE LIABLE FOR ANY INCIDENTAL OR CONSEQUENTIAL DAMAGES CAUSED BY ANY USE, DEFECT OR FAILURE OF THE PRODUCT, WHETHER THE CLAIM IS BASED ON WARRANTY, CONTRACT, TORT OR OTHERWISE.

The exclusions and limitations set out above are not intended to, and should not be construed so as to, contravene mandatory provisions of applicable law. If any part or term of this DISCLAIMER OF WARRANTY is held by any court of competent jurisdiction to be illegal, unenforceable or in conflict with applicable law, the validity of the remaining portion of the DISCLAIMER OF WARRANTY shall not be affected, and all rights and obligations shall be construed and enforced as if this DISCLAIMER OF WARRANTY did not contain the particular part or term held to be invalid.

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© 2006, 2007, 2011 Medtronic
FREESTYLE® and AOA® are registered trademarks of Medtronic, Inc.
Medtronic of Canada Ltd.